

crownworlddental.com | 866-497-3699

<b>R</b>			
Dr:	BruxZir Full-Strength	□ Full-Cast	NEW! Comfort3D Rite Splint
Account #:	(1,150 MPa)	□ Non-Precious	Bite Splint (3D-printed, hard)
	BruxZir Esthetic	Noble-Cast 45 (40% Au)	
Address:	(870 MPa)	(40 % Au) □ Noble-Cast 60	□ Comfort H/S Bite Splint (hard, with soft reline)
City:	🗌 IPS e.max (500 MPa)	(57.5% Au)	
	Obsidian All-Ceramic	🗌 White High Noble	
State	🗌 Obsidian PFM	(40% Au, 39.4% Pd)	
ZIP:	□ NP	White Noble Metal	
		(70.9% Pd)	
Phone:	□ WHN		
Dertie at ID (News et			
Due Date: Tooth #: Final Shade: Prep Shade:			
Due Date: Tooth #: Final Shade: Prep Shade:			
Patient ID/Name:	Doctor's signature require	ed:	
Due Date: Tooth #: Final Shade: Prep Shade:		ed: Date:	

## Your source for outstanding precision, consistency and value

- Lifetime warranty on BruxZir® Zirconia restorations
- **7-year warranty** on glass-ceramic, Obsidian<sup>®</sup> PFM, and full-cast restorations
- 6-month warranty on thermoformed appliances and splints

Crown World Dental Lab is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

For full warranty details, go to **crownworlddental.com/warranty**.

All rush cases must be prescheduled by calling **866-497-3699** before the case is shipped.

Time of pickup and delivery may affect turnaround time.

We accept all major credit cards.





**Signature Anterior Series** Choose the desired esthetic outcome for your patient SOFT-SQUARED SQUARE-TAPERED Exemplify boldness Show energetic and vouthfulness. professionalism. ☐TRAPEZOID Demonstrate Convey charm confidence. and softness. TRIANGLE-TAPERED

Rx must be enclosed with original case submission.